

DNA TOP-PCR Kit Protocol

Further information

- Technical assistance: www.topscibio.com

End Prep and adaptor ligation

- 1) Add the following components to a sterile nuclease-free PCR tube:

Fragmented DNA (2 ng)	6.6 μ L
End Modification Buffer	1.0 μ L
End Modification Enzyme	0.4 μ L
Total	8.0 μ L

- 2) Gently mix thoroughly and incubate at 20 °C for 30 min, following at 65 °C for 30 min and hold at 4°C.

- 3) Add the following components to the same PCR tube:

Adaptor	2 μ L
Ligation Mix (mix well before adding)	10 μ L
Total	20 μ L

- 4) Gently mix thoroughly and incubate at 20 °C for 4 hr and can be stored at -20 °C.

TOP-PCR reaction

- 5) Add the following components to a new sterile nuclease-free PCR tube:

Nuclease-free water	15 μ L
Ligation mixture from step 4 (mix well)	5 μ L
T3 oligo	5 μ L
2X PCR Mix	25 μ L
Total	50 μ L

[Note: the rest of ligation mixture can be stored at -20 °C as back-up]

- 6) Gently mix thoroughly and put in a PCR machine following the thermocycling condition:

Step	Temperature	Time	Cycles
Initial denaturation	98°C	30 sec	1
Denaturation	98°C	10 sec	15-30
Annealing	59°C	30 sec	
Elongation	72°C	1 min	
Final elongation	72°C	5 min	1
Hold	4°C	∞	1

- 7) Add 2 μ L of **Exonuclease** and gently mix thoroughly. Incubate at 37°C for 10 min, following at 80°C for 1 min and hold at 4°C.
- 8) Cleanup of PCR Reaction. Add 30 μ L of nuclease-free water or 0.1X TE buffer and store at -20 °C.